

## Weighing tables

- As the name indicates, the weighing tables are designed for equipment which are sensitive to vibration, such as scales. The large mass and damping elements made of rubber eliminate disturbances effectively.
- Materials are easy to clean.
- For Standard solutions we have selected high-quality and most commonly used options in terms of dimensions and materials.
- Custom solutions are always tailored to the customer's needs.



### Product information


Description	Product code and standard dimensions width-depth-height [mm]	Height [mm]	Weight approx. [kg]
Weighing table with electrically adjustable height	EVD-1200-550-S	Adjustment range from 650 to 940	170
Weighing table with electrically adjustable height	EVD-900-550-S	Adjustment range from 680 to 970	160
Weighing table with electrically adjustable height	EVD-550-550-S	Adjustment range from 680 to 970	100
Weighing table with a metal base	EVA-550-550-750	750	100
Weighing table with stone base	EVB-900-550-750	750	300
Weighing stone on a worktop	EVC-300-300-30	30	10

### Custom


In addition to Standard solutions, we also implement weighing tables according to customer's wishes.

Feel free to contact our experts and we will find suitable solutions together!

### VISU KALUSTE OY

 Eeliksentie 6, 34870 Visuvesi, Finland

 [visu@visu.fi](mailto:visu@visu.fi)

 +358 44 749 4940

 [visu.fi](http://visu.fi)




INNOSTAVIA YMPÄRISTÖJÄ  
TIETEELLE JA TUTKIMUKSELLE


## Features

<b>Installation</b>	Free-standing on the floor, must be placed away from other furniture and walls. With damping elements made of rubber.
<b>Materials</b>	Powder coated steel, acid-resistant solid laminate and polished granite.
<b>Features EVD (electrically adjustable)</b>	Load-bearing capacity 150 kg. Special wheels fully adjustable within the structure. When adjusted accordingly, the wheels functionally transform into table legs.
<b>By order</b>	A weighing table or its worktop can also be embedded into a laboratory desk or fume cupboard. Electrical outlets (3xSchuko IP21) can be added to the EVD electrically adjustable weighing tables.

### VISU KALUSTE OY

 Eeliksentie 6, 34870 Visuvesi, Finland

 visu@visu.fi

 +358 44 749 4940

 visu.fi