NORDIC BIOLABS PRESENTS THE

Christmas Calendar 2022

All offers in the Christmas Calendar are valid December 1st to 22nd .

Please think of the environment and order all your products at one occasion.

The offers can not be combined with other offers or discounts from Nordic Biolabs. Please refer to the Christmas Calendar for correct price when ordering. VAT and shipping cost will be added, but no small order fee or mosa will be charged.

Order tel. 32 52 88 00, e-mail info@nordicbiolabs.se

More information about the products can be found in our webshop - www.nordicbiolabs.dk.

1 dec

Vortex V-1 plus - the ultimate choice for mixing

40% discount



Two modes:

1. Continuous operation
2. Impulse operation (activated by pressing the cap with the tube's bottom)

V-1 Personal vortex
Art no: BS-010203-AAG

Now! 1380 DKK
List price 2301 DKK

biosan Medical–Biological Research & Technologies

2 dec

Nordic Biolabs Mini Centrifuge is compact, easy to use and comes with a host of extras

60% discount

- Small footprint and quiet operation
- Palm-shaped lid
- Change rotors with the push of a button

Maximum Speed: approx. 6,000 rpm
Maximum RCF: 2000 x g



Sprout Plus Mini-Centrifuge, Blue
Art no: HS120611

Now! 1512 DKK
List price 3780 DKK

Heathrow Scientific design | function | innovation

5 dec

FTA-1, Flask-trap aspirator, volume 1L

Designed for aspiration/removal of residual liquids (buffers, alcohols etc) from cell culture vessels, test tubes or microtubes, in cell culture or molecular biology (ex DNA/RNA precipitation).



40% discount

FTA-1, Flask-trap aspirator Volume 1L
Art no: BS-040108-AAG

Now! 2875 DKK
List price 4793 DKK

biosan Medical–Biological Research & Technologies

6 dec

PDS-250 is ready-to-use for eliminating DNA and RNA from suitable surfaces

35% discount

The use of PDS-250 both before and after PCR analysis is fast, easy and ideal to maintain a clean work area and thereby saves time and expenses.



PDS-250
Art no: BS-040107-DK

Now! 304 DKK
List price 468 DKK

biosan Medical–Biological Research & Technologies

nordic biolabs

www.nordicbiolabs.se info@nordicbiolabs.se ordertel 08-630 85 10 fax 08-630 09 05

7 dec

Restock with high quality water and buffer!
Your processes deserve it, not only during Christmas

40% discount




cytiva

CAPRICORN SCIENTIFIC

Art no.	Description	List price DKK	Now! DKK
Water For Injection, WFI			
SH31191.02	LM-HyPure WFI Quality Water	80	48
SH31191.03	LM-HyPure WFI Quality Water	115	69
Dulbecco's Phosphate Buffered Saline, D-PBS			
PBS-1A	D-PBS (1x) w/o Ca Mg w/o Phenol Red 500 mL	105	63
PBS-10XA	D-PBS (10x) w/o Ca Mg w/o Phenol Red 500 mL	205	123
PBS-2A	D-PBS (1x) w Ca Mg w/o Phenol Red 500 mL	102	61

8 dec

Sterile universal filter tips in rack



55% discount

Universal filter tips, RNase-, DNase-, Metal- and Pyrogen free. Good compatibility with major pipette brands like Eppendorf, Gilson etc.

- Sterilized by gamma irradiation SAL10-6 (ISO11137)
- Manufactured in a Class 100,000 cleanroom environment under ISO 13485:2016 and ISO 9001:2015 quality management system

deltalab

Art no.	Description	List price DKK	Now! DKK
D151020	20 µl Filter tip, in rack, 96x20 tips	2082	937
D151100	100 µl Filter tip, in rack, 96x20 tips	2082	937
D153200	200 µl Filter tip, in rack, 96x20 tips	2082	937
D151000	1000 µl Filter tip, in rack, 96x20 tips	2528	1138

9 dec

Oakton pH 550 Benchtop pH Meter Kit with Probe, Stand, and pH Buffers

Economical meter for benchtop pH measurement.

- Automatically recognizes calibration solutions
- Compact design saves valuable benchtop space
- Data logging allows you to store and recall up to 50 groups of data
- IP54-rated housing is dust- and spill-proof
- Bundled kit for convenience and savings

Oakton pH 550 Benchtop pH Meter Kit
with Probe, Stand, and pH Buffers
Art no: EW-35419-28

Now! 6470 DKK
List price 8091 DKK

20% discount



Cole-Parmer
an antylia scientific company

12 dec

The strongest pipette boy on the market



40% discount

- **Versatil usage:** fits all commercially available glass and plastic serological pipettes in the volume range 1-100 ml
- **Speed control:** an adjustment thumb wheel for setting the pipetting speed
- **Gravity dispensing:** dispensing drop-by-drop by applying a slight pressure on the dispense button

Pipettor Turbo Fix
Art no: 94700

Now! 2820 DKK
List price 4700 DKK

13 dec — Centrifuge tubes with benefits

- White writing surface on tube cap, works with any writing instrument
- Caps that cannot be over-tightened
- Graduations all the way down to 0.5 ml



- Continuous thread, external moulded graduations in 0,5 ml increments. Frosted area (55 x 10 mm).
- Blue cap in polyethylene with hermetical closure thanks to its internal elastic obturation
- Not autoclavable.
- Centrifugations force: 5000 xg



TPP

deltalab

Art no.	Description	List price DKK	Now! DKK
91050	Centrifuge tube, 50 ml, conical bottom, PP, 30x115 mm, 10x20 pcs	1479	890
429920-DL	Centrifuge tube, 15 ml, conical bottom, PP, sterile, 5x100 pcs	1036	700

14 dec — 3D "Sunflower Mini-Shaker" with adjustable three-dimensional smooth rotation of the platform

35% discount



The mini-shaker is designed for mixing blood samples, for minigel staining and destaining, sample washing, blot hybridization reactions.

3D, Sunflower Mini-Shaker
Art no: BS-010151-AAG

Now! 2632 DKK
List price 4049 DKK

biosan
Medical–Biological Research & Technologies

15 dec — Stands for pipettes to avoid contamination from benches and tables



Linear Stand for mechanical and electronical pipettes. The pipette stands hold your pipettes in an upright position while not in use to avoid contamination from benches and tables.

Charging Stand for Picus and eLINE 1 Place Pipette with Universal AC-Adaptor.

SARTORIUS

Art no.	Description	List price DKK	Now! DKK
725620	Linear stand for pipettes	845	590
730981	Charging stand for electronic pipette	2852	1800

16 dec — Portable spill kit that comes in a convenient dispenser box



The dispenser box contains 10 pads (41 x 51 cm). Designed for small spills of chemicals and hazardous liquids. Does not degrade or cause a dangerous chemical reaction with the absorbed liquid. Reusable lightweight plastic box that can be refilled once the sorbents have been used. Perfect to clean up chemical spills or for use in laboratories.

BRADY

Art no.	Size	List price DKK	Now! DKK
813910	10 pads/Kit	686	515
813842	100 pads	1706	1280

19 dec

Riplate PCR plate - affordable high quality

Riplate PCR plate i tested and works well on the widely used CFX96 DX System.

- 96 wells V-bottom, 150µl/well, skirted
- Material: PP (wells) / PC (frame)
- Standard format for use in automation
- Slightly raised well rims for better sealing
- Alphanumeric markings

40% discount



Riplate PCR plate 96, clear thin wall, low profile, 50 pcs/box
Art no: 43001-0305

Now! 1636 DKK
List price 2727 DKK

ritter MEDICAL

20 dec

Extract DNA in 3-8 min

The QuickExtract™ DNA Extraction Solution can be used to rapidly and efficiently extract PCR-ready genomic DNA from almost any sample type using a simple, one-tube protocol that takes only 3-8 minutes.

The QuickExtract method allows for the inexpensive processing of one to hundreds of samples simultaneously, without centrifugation, spin columns or the use of any toxic organic solvent. The method is also compatible with robotic automation.



QuickExtract DNA Extraction Solution 1.0
Art no: QE09050

Now! 3200 DKK
List price 5239 DKK

BIOSEARCH TECHNOLOGIES
Lucigen

21 dec

Spill control products for the lab



1. Kit for small spills of chemicals and hazardous liquids.
2. SpillFix granular is completely natural, producing no carcinogens or dust and has more than four times the absorption power of clay granular.
3. A chemical sorbent soc, filled with chemical meltblown polypropylene that absorbs chemicals and hazardous liquids.

BRADY

Art no.	Description	List price DKK	Now! DKK
813868	Chemical disposable spill kit	1079	809
149265	SpillFix Granular 3 Ltr. Jar	1040	780
813850	Soc, Ø 7.6 cm x 122 cm	1341	1006

22 dec

Reagent reservoirs, 25 and 50 ml



45% discount

- V-shaped well design – max reagent recovery
- Graduation marking with fill line on inner wall
- Drip free pour-off spouts in the corners for spill free pouring

Heathrow Scientific design | function | innovation

Art no.	Material	Sterile	Volume	Pack size	List price DKK	Now! DKK
HS120640	PS	Yes	25 ml	5x20 pcs	1945	1070
HS120642	PS	No	25 ml	100 pcs	783	430
HS120643	PS	Yes	50 ml	5x40 pcs	2287	1258
HS120645	PS	No	50 ml	100 pcs	865	475
HS120646	PVS	Yes	50 ml	5x20 pcs	1549	850

We wish you a Merry Christmas and a Happy New Year

The offers in the Christmas Calendar are valid until December 23. Please think of the environment and order all your products at one occasion. The offers can not be combined with other offers or agreements with Nordic Biolabs AB. All prices are ex VAT. Please refer to the Christmas Calendar in order for correct price when ordering. With reservation for misprints.

nordic biolabs

www.nordicbiolabs.se
Nordic Biolabs AB, Box 7293, 187 14 Täby,
Tlf. 32 52 88 00, Fax 32 52 88 79, info@nordicbiolabs.se